



Food Microbiology Proficiency Testing

ISO 17034 Accredited



spex.com

Phone: +1.732.549.7144 • +1.800.LAB.SPEX
Fax: +1.732.603.9647
spexsales@antylia.com



Spex and NSI Lab Solutions are
Antylia Scientific companies.

nsilabsolutions.com

Phone: +1.800.234.7837
Fax: +1.919.789.3019
nsi@nsilabsolutions.com

Food Microbiology Proficiency Testing

Innovative Controls to Support Data Quality

NSI Lab Solutions Food Microbiology Proficiency Testing (PT) studies are designed for your convenience and success.

- 4 studies yearly without order deadlines. You must choose a study number.
- PT Express™ quick turnaround PT studies available.
- Simple online reporting with no limits on the number of analysts or number of methods that you may report. Final reports returned within 21 days of the study close.

2022 Food Microbiology Study Schedule

Study Number	Open Date	Close Date
FM-0822	August 9	September 22
FM-1122	November 2	December 16



Certificate #: AT-1690



Certificate #: AR-1571



Certificate #: AP-1693-1



10002343 QM08
ISO 9001:2008

Food Microbiology Proficiency Testing Standards

Quantitative Indicators

A blended organism OT standard for quantitative determination of APC, total coliform, E. coli, S. aureus, enterobacteriaceae, and yeast/mold in the range of 100–250,000 CFU/g after hydration. Supplied with hydration fluid.

Type	Part #
Buffer Matrix (2 pack)	FMPT-001
Buffer Matrix (5 pack)	FMPT-001-5
Buffer Matrix (10 pack)	FMPT-001-10
Meat Matrix (2 pack)	FMPT-001M
Meat Matrix (5 pack)	FMPT-001M-5
Meat Matrix (10 pack)	FMPT-001M-10
Dairy Matrix (2 pack)	FMPT-001D
Dairy Matrix (5 pack)	FMPT-001D-5
Dairy Matrix (10 pack)	FMPT-001D-10

Also available as PT Express samples.

Quantitative Lactic Acid Bacteria

A pure culture conveniently supplied in duplicate for quantitative determination of lactic acid producing bacteria, in the range of 100–250,000 CFU/g after hydration. Supplied with hydration fluid.

Type	Part #
Buffer Matrix	FMPT-004
Meat Matrix	FMPT-004M
Dairy Matrix	FMPT-004D

Also available as PT Express samples.

Quantitative Bacillus cereus

A pure culture conveniently supplied in duplicate for quantitative determination of B. cereus in the range of 100–250,000 CFU/g after hydration. Supplied with hydration fluid.

Type	Part #
Buffer Matrix	FMPT-003
Meat Matrix	FMPT-003M
Dairy Matrix	FMPT-003D

Also available as PT Express samples.

Quantitative Campylobacter

A pure culture conveniently supplied in duplicate for quantitative determination of campylobacter, in the range of 100–250,000 CFU/g after hydration. Supplied with hydration fluid.

Type	Part #
Buffer Matrix	FMPT-016
Meat Matrix	FMPT-016M
Dairy Matrix	FMPT-016D

Also available as PT Express samples.

2022 Food Microbiology Study Schedule

Study Number	Open Date	Close Date
FM-0822	August 9	September 22
FM-1122	November 2	December 16

Food Microbiology Proficiency Testing Standards (continued)

Quantitative Psuedomonas

A pure culture of *P. aeruginosa* conveniently supplied in duplicate for quantitative determination of psuedomonas, in the range of 100–250,000 CFU/g after hydration. Supplied with hydration fluid.

Type	Part #
Buffer Matrix	FMPT-015
Meat Matrix	FMPT-015M
Dairy Matrix	FMPT-015D

Also available as PT Express samples.

Qualitative Listeria monocytogenes Set

A three sample set for qualitative identification of *L. monocytogenes*. Set contains *L. monocytogenes* plus two non-pathogen, non-listeria organisms. Report present or absent for *L. monocytogenes* for each sample in the set. Three different series available per study. Supplied with hydration fluid.

Type	Part #
Water Matrix	FMPT-006
Meat Matrix	FMPT-006M
Dairy Matrix	FMPT-006D

Also available as PT Express samples.

Qualitative Salmonella Set

A three sample set for qualitative identification of salmonella, spp. Set contains salmonella and non-pathogenic, non-salmonella organisms. Report present or absent for salmonella spp. for each sample. Three different series available per study. Supplied with hydration fluid.

Type	Part #
Water Matrix	FMPT-008
Meat Matrix	FMPT-008M
Dairy Matrix	FMPT-008D

Also available as PT Express samples.

Qualitative STEC

A single sample conveniently supplied in duplicate containing at least 1 of 6 STECs for identification. Supplied with hydration fluid.

Type	Part #
Water Matrix	FMPT-009
Meat Matrix	FMPT-009M
Dairy Matrix	FMPT-009D

Also available as PT Express samples.

2022 Food Microbiology Study Schedule

Study Number	Open Date	Close Date
FM-0822	August 9	September 22
FM-1122	November 2	December 16

Food Microbiology Proficiency Testing Standards (continued)

Quantitative Clostridium perfringens

A pure culture of C. perfringens conveniently supplied in duplicate for quantitative determination of C. perfringens. Supplied with hydration fluid.

Type	Part #
Buffer Matrix	FMPT-014
Meat Matrix	FMPT-014M
Dairy Matrix	FMPT-014D

Also available as PT Express samples.

Qualitative Pathogens Set

A unique three sample set for qualitative identification of S. enterica, L. monocytogenes and E. coli O157:H7. Report present or absent for each organism of interest. Three different series available per study. Supplied with hydration fluid.

Type	Part #
Water Matrix	FMPT-005
Meat Matrix	FMPT-005M
Dairy Matrix	FMPT-005D

Also available as PT Express samples.

Qualitative Listeria Set (Non-Pathogenic)

A three sample set for qualitative identification of listeria spp. Set contains non-pathogenic listeria strain and non-pathogenic, non-listeria organisms. Report present or absent for listeria spp. for each sample in the set. Three different series available per study. Supplied with hydration fluid.

Type	Part #
Water Matrix	FMPT-007
Meat Matrix	FMPT-007M
Dairy Matrix	FMPT-007D

Also available as PT Express samples.

Qualitative Campylobacter

A three sample set for qualitative identification of campylobacter spp. Set contains C. jejuni or coli and non-pathogen/non-campylobacter facultative anaerobes. Report present or absent for campylobacter spp for each sample. Supplied with hydration fluid.

Type	Part #
Water Matrix	FMPT-013
Meat Matrix	FMPT-013M
Dairy Matrix	FMPT-013D

Also available as PT Express samples.

2022 Food Microbiology Study Schedule

Study Number	Open Date	Close Date
FM-0822	August 9	September 22
FM-1122	November 2	December 16

Food Microbiology Proficiency Testing Standards (continued)

Qualitative *Listeria monocytogenes* Set

A three sample set for qualitative identification of *L. monocytogenes*. Set is designed to mimic typical FDA BAM and USDA MLG sample handling procedures. Each set contains 3 x 25 grams of sterile matrix (beef powder or skim milk powder) and three individually packaged lyophilized microorganism pellets. At least one of three of the sample pellets will be positive for *L. monocytogenes*.

To use, transfer 25 grams of sterile matrix to stomacher bag, add 225 mL of laboratory supplied enrichment broth and then add a single sample pellet to the stomacher bag. Process and analyze according to your normal laboratory procedures. Report present or absent for each sample set. Three different series available per study.

Type	Part #
Meat Matrix	FMPT-022M
Dairy Matrix	FMPT-022D

Also available as PT Express samples.

Qualitative STEC Set

A three sample set for qualitative identification of STEC. Set is designed to mimic typical FDA BAM and USDA MLG sample handling procedures. The PT set will contain one of the following: O26, O45, O103, O111, O121 or O145. The PT set can be used to report Total STEC or the specific serotype. Each set contain 3 x 25 grams sterile matrix (beef powder and skim milk powder) and three individually packaged lyophilized microorganism pellets. At least one of the three samples will be positive for STEC.

To use, transfer on 25 gram sterile matrix to a stomacher bag, add 225 mL of your in-house enrichment broth and then add a single sample pellet to the stomacher bag. Process and analyze according to your normal laboratory procedures. Report present or absent for each sample set. Three different series available per study.

Type	Part #
Meat Matrix	FMPT-024M
Dairy Matrix	FMPT-024D

Also available as PT Express samples.

2022 Food Microbiology Study Schedule

Study Number	Open Date	Close Date
FM-0822	August 9	September 22
FM-1122	November 2	December 16

Food Microbiology Proficiency Testing Standards (continued)

Qualitative Salmonella Set

A three sample set for qualitative identification of Salmonella spp. Set is designed to mimic typical FDA BAM and USDA MLG sample handling procedures. Each set contains 3 x 25 grams of sterile matrix (beef powder or skim milk powder) and three individually packaged lyophilized microorganism pellets. At least one of three of the sample pellets will be positive for Salmonella.

To use, transfer on 25 gram sterile matrix to a stomacher bag, add 225 mL of your in-house enrichment broth and then add a single sample pellet to the stomacher bag. Process and analyze according to your normal laboratory procedures. Report present or absent for each sample set. Three different series available per study.

Type	Part #
Meat Matrix	FMPT-023M
Dairy Matrix	FMPT-023D

Also available as PT Express samples.

Qualitative E. coli O157:H7 Set

A three sample set for qualitative identification of E. coli O157:H7. Set is designed to mimic typical FDA BAM and USDA MLG sample handling procedures. Each set contains 3 x 25 grams of sterile matrix (beef powder or skim milk powder) and three individually packaged lyophilized microorganism pellets. At least one of the three pellets will be positive for E. coli O157:H7.

To use, transfer on 25 gram sterile matrix to a stomacher bag, add 225 mL of your in-house enrichment broth and then add a single sample pellet to the stomacher bag. Process and analyze according to your normal laboratory procedures. Report present or absent for each sample set. Three different series available per study.

Type	Part #
Meat Matrix	FMPT-025M
Dairy Matrix	FMPT-025D

Also available as PT Express samples.

Environmental Swab - Listeria spp.

A 5 sample panel for qualitative listeria identification of Listeria spp. in environmental swabs. Each panel is supplied with 5 inoculated swabs and five 4"x4" sterile swabbing surfaces. At least 2 of 5 inoculated swabs will be positive for Listeria spp. Acceptable evaluation is 4 of 5 correct with no false negatives. Listeria species utilized for this panel is Listeria ivanovii.

Part #
FMPT-018

Also available as PT Express samples.

2022 Food Microbiology Study Schedule

Study Number	Open Date	Close Date
FM-0822	August 9	September 22
FM-1122	November 2	December 16

Food Microbiology Proficiency Testing Standards (continued)

Environmental Swab - E. coli O157:H7

A five sample panel for qualitative identification of E. coli O157:H7 in environmental swabs. Each panel is supplied with 5 inoculated swabs and five 4"x4" sterile swabbing surfaces. At least 2 of 5 inoculated swabs will be positive for E. coli O157:H7. Acceptable evaluation is 4 of 5 correct with no false negatives.

Part #
FMPT-021

Also available as PT Express samples.

Environmental Swab - Salmonella spp.

A five sample panel for qualitative identification of Salmonella spp. in environmental swabs. Each panel is supplied with 5 inoculated swabs and five 4"x4" sterile swabbing surfaces. At least 2 of 5 inoculated swabs will be positive for Salmonella spp. Acceptable evaluation is 4 of 5 correct with no false negatives.

Part #
FMPT-019

Also available as PT Express samples.

Environmental Swab - Listeria monocytogenes

A five sample panel for qualitative identification of L. monocytogenes in environmental swabs. Each panel is supplied with 5 inoculated swabs and five 4"x4" sterile swabbing surfaces. At least 2 of 5 inoculated swabs will be positive for L. monocytogenes. Acceptable evaluation is 4 of 5 correct with no false negatives.

Part #
FMPT-020

Also available as PT Express samples.

Environmental Swab - Listeria sp.

A five sample panel for qualitative identification of Listeria spp. in environmental swabs. Each panel is supplied with 5 inoculated swabs and five 4"x4" sterile swabbing surfaces. At least 2 of 5 inoculated swabs will be positive for Listeria spp. Acceptable evaluation is 4 of 5 correct with no false negatives. Listeria species utilized for this panel is Listeria ivanovii.

Part #
FMPT-018

Also available as PT Express samples.

2022 Food Microbiology Study Schedule

Study Number	Open Date	Close Date
FM-0822	August 9	September 22
FM-1122	November 2	December 16

spex.com

Phone: +1.732.549.7144 • +1.800.LAB.SPEX
 Fax: +1.732.603.9647
 spexsales@antylia.com



Spex and NSI Lab Solutions are Antylia Scientific companies.

nsilabsolutions.com

Phone: +1.800.234.7837
 Fax: +1.919.789.3019
 nsi@nsilabsolutions.com