

Food Chemistry Proficiency Testing

ISO 17043 Accredited

Innovative Controls to Support Data Quality

NSI Lab Solutions Food Chemistry Proficiency Testing (PT) studies are designed for your convenience and success.

- 4 studies yearly without order deadlines. You must choose a study number.
- PT Express™ quick turnaround PT studies available.
- Simple online reporting with no limits on the number of analysts or number of methods that you may report. Final reports returned within 21 days of the study close.



Qualitative Allergens Panels

Each panel includes 3 samples with at least 1 sample containing the allergen of interest at a level close to regulatory threshold. Verified to work with various test technologies. Each vial contains approximately 10 grams of material, 3 distinct series of each panel are available each study.

Type	Part #
Qualitative Gluten	FCPT-007
Qualitative Peanut	FCPT-008
Qualitative Egg	FCPT-009
Qualitative Milk	FCPT-010
Qualitative Crustacean	FCPT-011
Qualitative Soy	FCPT-012

Also available as PT Express samples.

pH and Titratable Acidity in Dairy COMING SOON

The PT material is typically a skim milk. Analyze for pH and titratable acidity. Supplied in 2 x 100 mL bottles.

Part #
FCPT-013

Qualitative Gluten on Surface COMING SOON

A five sample panel for qualitative identification of Gluten on Surface. Each panel consists five 4"x4" plates that have been pre-inoculated. At least two of the five inoculated plates will be positive for Gluten on Surface. Acceptable evaluation is 4 of 5 correct with no false negatives. This product is currently not covered under our ANAB scope.

Part #
FCPT-016

Gluten in Fermented Product COMING SOON

A quantitative single sample of lyophilized soy sauce or Guinness Stout. Applicable for r-Biopharm 7021 or equivalent methods. This product is currently not covered under our ANAB scope.

Part #
FCPT-018

Gluten in Food Product COMING SOON

A quantitative single sample of gluten in rice flour in the range of 10-200 mg/kg. Applicable for Neogen Veratox, 3M and r-Biopharm methods. This product is currently not covered under our ANAB scope.

Part #
FCPT-021

Pesticides in Lettuce COMING SOON

A lyophilized 10 gram slurry, fortified with pesticides of interest. Supplied in duplicate. This product is currently not covered under our ANAB scope.

Part #
FCPT-014

Pesticides in Grapes COMING SOON

A lyophilized 10 gram slurry, fortified with pesticides of interest. Supplied in duplicate. This product is currently not covered under our ANAB scope.

Part #
FCPT-015

10002343 QM08



NSI Lab Solutions

ISO 9001:2015



Mycotoxins in Corn Meal
COMING SOON

A single quantitative sample for determination of various mycotoxins.

Part #
FCPT-017

Heavy metals in Baby Food-Sweet Potato
COMING SOON

A wet sweet potato fortified with heavy metals of interest at levels 2-5 times action levels.

Part #
FCPT-020

Heavy Metals in Baby Cereal
COMING SOON

A rice cereal fortified with heavy metals of interest at levels 2-5 time action levels.

Part #
FCPT-019

2023 Food Chemistry Study Schedule

Study Number	Open Date	Close Date
FM-0223	Feb. 15, 2023	March 31, 2023
FM-0523	May 10, 2023	June 23, 2023
FM-0823	Aug. 9, 2023	Sept. 22, 2023
FM-1123	Nov. 7, 2023	Dec. 21, 2023

10002343 QM08



NSI Lab Solutions

ISO 9001:2015

