




## Mycotoxins

- Mycotoxins are ISO 17034 accredited
- Mix includes 5 of the most commonly tested mycotoxins
- Custom standards are available



**spex.com**

Phone: +1.732.549.7144 • +1.800.LAB.SPEX  
Fax: +1.732.603.9647  
spexsales@antylia.com

Connect with us



Spex CertiPrep is an  
Antylia Scientific company.  
Find out more at [antylia.com](http://antylia.com).



# Mycotoxins

Mold and fungi are all over the world and are of major concern in agricultural products. Fungi and their spores can spread throughout crops and food stores tainting crops. Fungi can proliferate into colonies with the right environmental conditions which destroy crops by consuming them or poison them with toxic secondary metabolites called mycotoxins.

Spex CertiPrep offers an analytical standard mix of five of the most commonly tested mycotoxins (four aflatoxins and ochratoxin A) for your mycotoxin testing needs. Our mycotoxins are ISO 17034 accredited. Custom standards are also available. Contact us at +1.732.549.7144 or via email at [spexsales@antylia.com](mailto:spexsales@antylia.com).

Mycotoxins Standard Mix					
Compound	CAS #	Concentration	Matrix	Volume	Part #
Aflatoxin B1	1162-65-8	10 µg/mL of each compound	Acetonitrile	1 mL	MYCO-1
Aflatoxin B2	7220-81-7				
Aflatoxin G1	1165-39-5				
Aflatoxin G2	7241-98-7				
Ochratoxin A	303-47-9				



Organic Certified  
Reference Materials



Supplied with a  
Certificate of Analysis



ISO Accredited  
Standards

**spex.com**

Phone: +1.732.549.7144 • +1.800.LAB.SPEX  
Fax: +1.732.603.9647  
[spexsales@antylia.com](mailto:spexsales@antylia.com)

4772W

Connect with us



Spex CertiPrep is an  
Antylia Scientific company.  
Find out more at [antylia.com](http://antylia.com).

